



2° HEEL WEDGES

MATERIALS

- Durable, dense EVA foam, Blue
- Adhesive on bottom with removable backing

FEATURES & BENEFITS

- Provides 2 degrees of medial posting for greater correction
- Helps prevent over-pronation and flattening of the foot
- Adhesive bottom to help hold heel wedges in place
- Provides enhanced support for those who need it
- Sold as pair

IDEAL FOR: Over-pronation, flat feet, fallen arches, arch and heel pain

Part No.	Size	US Shoe Size	UPC Code
5095-01S	Small	M 3-8, W 5-10	749852100185
5095-01L	Regular	M 9-14, W 11-16	749852100284

APPLICATION INSTRUCTIONS

For use with Powerstep® Insoles: Remove adhesive backing, noting the Right (R) and Left (L) markings. Apply firmly to the underside of the Powerstep® orthotic heel platform with the curved edge aligned with the curve of the heel. Adjust position as necessary for comfort.

Not intended for persons with severe foot deformities or circulatory disorders resulting in infections or foot ulcerations. Always consult a medical professional when experiencing persistent pain, discomfort or swelling.



888.237.3668 · 513.825.1888

Email: info@powersteps.com

www.powerstep.com

Stable Step, LLC, West Chester, OH

Made in the USA with US and imported materials.

Powerstep® is a registered trademark of Stable Step, LLC.